

Penn State New Kensington
Biomedical Engineering Technology

Electrical Safety

Myron D. Hartman MS, SASHE, CCE, CRES, CBET

Outcome 1, 2 & 6

BET 204W – Electrical Safety

Classes of Medical Devices

- Class G - General Areas
- Class H - Critical Care Areas
- Class W - Wet Locations

BET 204W – Electrical Safety

Electrical Safety

Leakage current – Where does it come from?

Impedance (Z) consists of three components

X_L – Inductance Reactance

X_C – Capacitance Reactance

R – Resistance

page 715

Carr-Brown

$$Z = \sqrt{R^2 + (X_L - X_C)^2}$$

BET 204W – Electrical Safety

Electrical Safety

Leakage current

Micro & Macro Shock

Micro Shock or Cardiac Shock - is defined as a low value current (ua), which passes directly through the heart via a needle or catheter in an artery or a vein

- Less than 1ma of current.
- Current can still have effect on patients with invasive equipment, such as: pacemakers, swan-ganz catheters, invasive blood pressure lines, cardiac cath lines, etc.

page 714

Carr-Brown

BET 204W – Electrical Safety

Electrical Safety

Leakage current

Micro & Macro Shock

Macro Shock – is defined as a high-value current (ma), which passes arm to arm through the body by (skin) contact with a voltage source. There must be two points of body contact. The resulting current eventually passes through the heart and may cause ventricular fibrillation of death.

- Body can sense shock.
- 1ma - Threshold of Perception
 - 10-20ma - Let go Current
 - 50ma - Pain
 - 100-300ma - Ventricular Fibrillation
 - 6A - Total Myocardial Contraction
 - >6A - Burns and Respiratory Paralysis

page 714

Carr-Brown

BET 204W – Electrical Safety

Electrical Safety

Leakage current

High frequency currents – The primary concern of electrical shock is with 60 hz currents. Medical equipment today may produce leakage current at higher frequencies, due to harmonics, carrier frequencies and switching mode power supplies.

Frequencies above 100 khz have current limits increased 100 times that of low frequency currents

Frequencies from 1 khz to 100 khz have their limits increased on a linear scale up to 100 times that of the low frequencies.

As the frequency is increased, there is less effect of shock and increased effect of RF burns

BET 204W – Electrical Safety


Electrical Safety

Leakage current

AAMI filter – Due to multiple frequencies that could exist in leakage current the AAMI filter is used to pass only the low frequencies of interest. High frequencies are shunted to ground, This filter provides a weighted risk current.

page 628

PBETMI.





BET 204W – Electrical Safety

Electrical Safety

Leakage current

Standards

- AAMI – Association for the Advancement of Medical Instrumentation
- NFPA 99 – National Fire Protection Agency – Health Care Facilities (2005 Edition)



BET 204W – Electrical Safety

Electrical Safety

Leakage current

Current limits – are the safe levels of currents produced by the medical equipment, including all cables.

Chassis current – current on the metal cabinet

Isolated Patient Current – current on any patient connection that has an isolated circuit

Non-Isolated Patient Current – current on any patient connections that is not isolated (older equipment)

BET 204W – Electrical Safety

Electrical Safety

Leakage current

Current limits – Below are the **AAMI** recommended leakage current limits.

- Chassis, Non-Patient Contact - 500 ua
- Chassis, Patient Contact - 100 ua
- Patient Leadwire, Isolated - 10 ua
- Patient Leadwire, Sink Current - 20 ua
- Patient Leadwire, Non-Isolated - 50 ua

page 96

1985 AAMI Std.

BET 204W – Electrical Safety

Electrical Safety

Leakage current

Current limits – Below are the **NFPA** recommended leakage current limits.

- Chassis, Fixed Equipment – 5 ma
- Chassis, Patient Contact - 300 ua
- Patient Leadwire (all to G), Isolated - 100 ua
- Patient Leadwire (single to G, w/G), Isolated - 10 ua
- Patient Leadwire (single to G, wo/G), Isolated - 100 ua
- Patient Leadwire, Sink Current - 50 ua
- Patient Leadwire, Non-Isolated - 50 ua

page 632

PBETMI.

BET 204W – Electrical Safety

Electrical Safety

Leakage current

Current limits - are measured in the following configurations:

- | | | |
|---------------------------------|--------|---------|
| • Device ground | OPEN | CLOSED |
| • Device power | ON | OFF |
| • Device polarity (hot/neutral) | NORMAL | REVERSE |

Record the highest measured leakage current and the configuration (ON/Off, Ground OPEN/CLOSED and Polarity NOR/REV)

BET 204W – Electrical Safety

Electrical Safety

Leakage current

120 Volt Isolation test (sink current) – This is a test performed on the isolation of circuits of medical equipment, to verify that an external current can not pass from the patient back into the medical equipment.

A potential of 120 Volts rms, 60 hz, is applied through a 120 K ohm resistance, to each patient connection of the medical device. The sink current that is measured is the amount that passes back into the medical device.

The test is performed with the medical device powered on and off.

page 633

PBETMI.

BET 204W – Electrical Safety

Electrical Safety

Leakage current

120 Volt Isolation test (sink current)

page 98
1985 AAMI Stds.

page 633
PBETMI.




Figure 6. Patient sink current test circuit.

Note: The 120-kilohm resistance is intended to protect the test operator.

BET 204W – Electrical Safety

Electrical Safety

Grounding

Why it is necessary – Provides the least amount of resistance for leakage current to flow back to the power source.

Current will always take the least path of resistance.

Leakage current exists in all electrically powered devices, therefore a ground is required to remove the current.

Ground conductor is also required for faults or failures with the device so currents have a safe path to follow.

BET 204W – Electrical Safety

Electrical Safety

Grounding

Where it is measured – from the device chassis to the ground pin of the power plug.

If possible, the plug should be removed from the receptacle with the power off.

If the device is “hard wired,” the ground can still be measured from the metal chassis to a ground on the electrical distribution system (power off).

BET 204W – Electrical Safety

Electrical Safety

Grounding

Measurement limits – the maximum grounding resistance should not be greater than 0.5 ohms.

While measuring, the power cord should be twisted and moved at the strain relief to the device and at the power plug to check for loose connections.

Redundant grounds should be removed if the device can not be unplugged.

BET 204W – Electrical Safety

Power Plugs & Outlets

- Hospital Grade - Green Dot
- Cut Wires for hot/neutral shorter
- Ground Pin Longer - first to make, last to break
- Wire Colors:
 - Ground – Green (US) or Green-Yellow Stripe (Europe)
 - Hot – Black (US) or Brown (Europe)
 - Neutral – White (US) or Blue (Europe)



BET 204W – Electrical Safety

Wiring Devices

120 Volt Plugs 3 prong (15 and 20 amp)

Installation

Hot and Neutral wires cut shorter (1/4 inch)



BET 204W – Electrical Safety

Power Plugs & Outlets

- Normal Power Outlets - Ivory/**Brown**
- Normal/Emergency Outlets - **Red**
- Isolated Ground Outlets - **Orange**
- Surge Protected – **Blue**
- Tamper Proof - Pediatric & Psychiatric
- Ground Fault Circuit Interrupters (GFCIs) – Class W
- Cover Plates & Circuit Identification
- Hospital Grade Rating

pages 763-764

PBETMI.

Power Plugs & Outlets

- Ground Fault Circuit Interrupter (GFCI)
Wet Locations, Shock Hazard Elimination, No
Life Support Equipment, Trip Point 5-6 ma.

- Line Isolation Monitors (LIMS)

Shock Hazard Elimination, Trip Point 5-20 ma,
Isolation Transformer

page 756

PBETMI.

Isolated Power and Line Isolation Monitors (LIMs)

Why they are used

- For added protection to prevent electrical shock to the patient and the users of medical equipment
- Patients that have invasive procedures or catheters connected to their heart are at a greater risk of electrical shock
- Isolated power systems are simply “isolated ungrounded electrical distribution systems”
- They are not required by code for new installations, but any existing must be maintained, tested and calibrated


page 760

PBETMI.

BET 204W – Electrical Safety

Isolated Power and Line Isolation Monitors (LIMs)


Photos of Systems



BET 204W – Electrical Safety

Isolated Power and Line Isolation Monitors (LIMs)


Schematics



Page 58

Bio-Tek

Electrical Safety Guide



BET 204W – Electrical Safety

Isolated Power and Line Isolation Monitors (LIMs)

Transformer ratings

Typical transformer ratings for Isolated Power Systems are rated in KVA, with the most common being 5 and 7.5 KVA


The transformer turns ratio is simply a 1:1 with 120 V input and 120 V secondary

Most operate at 208 V input and 120 V secondary, which is a 1.73:1 turns ratio

Power in must equal power out, so the KVA rating can be used to calculate current on the secondary, ie: 7.5 KVA @ 120 V single phase will provide 62.5 amps of current ($7500/120$)

BET 204W – Electrical Safety

Isolated Power and Line Isolation Monitors (LIMs) Current loops



Page 59

Bio-Tek

PENNSTATE

Electrical Safety Guide

BIOMEDICAL
Engineering Technology

BET 204W – Electrical Safety

Isolated Power and Line Isolation Monitors (LIMs)

Testing

Secondary voltage – measured from X1 to X2

120 V RMS is the expected, $\pm 5\%$ (114 – 126 V)

60 Hz is the expected, $\pm 0.5\%$ (59.7 – 60.3 Hz)

Macro current A & B

Simulated short circuit from X1 or X2 to ground and measuring the maximum current that would flow in the system

LIM trip point

Set point should be about 5-6 ma

Checked with variable resistor from X1 or X2 to ground

LIM accuracy

Verifying accuracy of gauge or digital display against calibrated analyzer


Alarms

BET 204W – Electrical Safety

Isolated Power and Line Isolation Monitors (LIMs) Testing

Macro current A & B

Simulated short circuit from X1 or X2 to ground and measuring the maximum current that would flow in the system




BET 204W – Electrical Safety

Isolated Power and Line Isolation Monitors (LIMs)

Testing

LIM trip point



BET 204W – Electrical Safety


Isolated Power and Line Isolation Monitors (LIMs)

Outlet wiring

X1 –Orange

X2 –Brown

Ground – Green



BET 204W – Electrical Safety


Portable Isolation Transformers

Used for multiple systems in single configuration

Surgical video systems

Portable video equipment

Computers, printers, video displays, network devices,...



Portable Isolation Transformers

Used to reduce leakage current levels

The isolation transformer creates a secondary voltage that is not referenced to ground

Current will not flow from the secondary winding to ground, therefore leakage currents are eliminated

Typical sizes

0.25 KVA – 5.0 KVA

Transformers can be purchased with systems or added on as needed

BET 204W – Electrical Safety

120 Volt Receptacle Testing

Voltage & Frequency

Measured from Neutral to Hot on receptacle

120 V RMS is the expected, $\pm 5\%$ (114 – 126 V)

60 Hz is the expected, $\pm 0.5\%$ (59.7 – 60.3 Hz)

Wiring Polarity

The Hot connection is connected to the hot on the electrical distribution system

The Neutral connection is connected to the neutral on the electrical distribution system, which is connected to the ground

Pin Tension

At least 4 ounces of pull on each blade

Ground to Ground resistance

Must be less than 0.1 ohms

Ground to Neutral resistance

Must be less than 0.1 ohms

Faceplate cracked/broken

page 631

PBETMI.

Electrical Distribution Systems

Hospitals have 2 primary feed circuits that must take the full load of the building

These are connected by tie breaker that works automatically or manual

Under normal conditions, $\frac{1}{2}$ of the building is on one circuit and the other $\frac{1}{2}$ is on the other circuit

page 754

PBETMI.

Emergency Generator Testing

Emergency generator must supply power if utility power is lost or if there is an internal distribution failure. Power must be supplied to: ICCU, CCU, Ors, ED, Neuro, L/D, and all critical care areas

Testing - Emergency generators must be tested by starting them once a week and then placing them under load once a month

page 755

PBETMI.

Emergency Generator Testing

Power must be supplied within 10 seconds from time of utility interruption



page 755
PBETMI.

BET 204W – Electrical Safety

Electrical Distribution Systems

Automatic Transfer Switches (ATS) switch critical loads from either utility or emergency power.

Open Transition – Switch opens from one source (utility) then closes other source (emergency). This creates an interruption of power for a fraction of a second, but long enough to cause computers to reset.

Closed Transition – The switch parallels the power from the utility with the emergency generator for several m sec, then opens one source. No interruption of power occurs to loads down stream

page 753-754

PBETMI.

BET 204W – Electrical Safety

Electrical Distribution Systems

Automatic Transfer Switches (ATS)

